



## **Methods for Detection of Microorganisms False and False (Specificity and sensitivity)**

**Joan B. Rose**

# Methods and Microbes

- No method for concentration and detection is 100% efficient
- Microbes are individual particulates not solutes, and are not necessarily evenly distributed in a given media (air,water,soil)
- Concentrate
- Purify
- Separate
- Detect & Quantify (culture, microscopic, indirect)
- Determine viability (CULTURE)
- Determine Hazard (carrying virulence genes)

# Issues in Detection

- Sensitivity – Detection Volume and presence of interfering substances that reduce assay sensitivity
- Specificity – is it the right organism or group of organisms
- Quantification – precision
- TWO ways to test new approaches
  - Compare to golden standard
  - Seeded recovery methods (ratio of recovered concentration to a known concentration added to the matrix)

# Methods and Microbes

- False positive: A test result that indicates a positive status when the true status is negative
- False negative: A test result that indicates a negative status when the true status is positive
- A variety of probability estimates can be computed from the information organized in a two-way table given in Table 1.
- *Sensitivity*: Probability of a positive test result (or presence of the symptom) given the presence of the disease
  - *Specificity*: Probability of a negative test result (or absence of the symptom) given the absence of the disease
  -
- *Predictive value positive*: Probability that a subject has the disease given that a subject has a positive screening test result (or has the symptom)
  -
- *Predictive value negative*: Probability that a subject does not have the disease, given that the subject has a negative screening test result (or does not have the symptom)


# Cultural/Standard Methods for the Detection of Pathogens in the Environment

- Virus
  - Concentration/separation
  - Cell culture (one cell line does not detect all)
  - Serology (identification)
- Bacteria
  - Concentration/separation
  - Enrichment media
  - Selective media
  - Biochemical tests, serology, immunochemical

# Methods for the Detection of Pathogens

- Protozoa
  - Concentration or elution
  - Purification
    - Differential centrifugation
    - Immune Magnetic Separation
  - Stain with monoclonal antibodies
  - Observe under UV light
  - Examine for characteristics
    - Shape and size
    - Internal structures

# Detection Targets



# Animal Virus Structure




Figure 2



# Typical Sample Volumes for Water

- Bacteria
  - 100 ml
- Viruses
  - Raw Sewage 1-5 liters
  - Treated sewage 40 liters
  - Surface waters 400 liters
  - Drinking water 400 to 2000 liters
- Protozoa
  - Treated sewage 4 liters
  - Surface waters 10 to 100 liters
  - Drinking water 10 to 100 liters

# Time and Cost for Assay of Enteric Pathogens in Water

- Bacteria
  - 2 to 5 days; \$40 to \$200
- Virus
  - 14 to 60 days; \$500 to \$1,200
- Protozoa
  - 2 to 3 days; \$250 to \$450

# Enumeration of Microbes

- Quantitative measurements
- Quantal measurements (presence/absence)
- CFU          Colony Forming Unit
- PFU          Plaque Forming Unit
- MPN          Most Probable Number (statistical method to estimate concentrations with multiple assays at dilutions to extinction [0])




PHONE

## Steps Involved in the Detection of Waterborne Enteric Viruses and Parasites




# Concentration of viruses from water







# Noninfected monkey kidney cells





# Monkey cells infected with poliovirus





# Principles of Protozoan Detection in the Environment

- Filtration of large volumes of water
- Concentration
- Purification
- Method of detection coupled with an ability to quantify the desired microorganisms

# Parasite Concentration

- Pall Gelman Envirocheck Sampling Capsule
- Used for US EPA method 1622
  - Protozoa detection:  
*Cryptosporidium*  
and *Giardia*




Envirocheck Sampling Capsule

# Methodology

## (filtration/concentration with Filta-Max)

- After a volume of water is passed through a filter cartridge, the eluted solution is then concentrated



Filta-Max™


# Methodology

- The concentrate is then purified
- This is accomplished by immunomagnetic separation (IMS)




**Puri-Max™**

# Purification by IMS



# Detection by IF

- Detection method is based on an indirect immunofluorescent antibody (IFA) stain





## Examination



Read slide under UV epifluorescent microscope  
Presumptive



*Cryptosporidium*  
4-6 mm



*Giardia*  
8-12 mm

Confirm  
using Differential Interference Contrast



see up to four sporozoites



see nuclei,  
axoneme &  
median bodies


# Cell Culture-*Cryptosporidium*




Take purified sample  
expose to 10% bleach for  
10 min to inactivate  
viruses and bacteria,  
algae, fungi



Wash and perform  
excystation





Inoculate sample onto  
HCT-8 cells



Examine cells for  
evidence of infection and  
growth of *Cryptosporidium*  
life stages

- microscope (bright field/IFA)
- ELISA
- PCR












# Aerosol Sampling

- Impingers
- Impact- Anderson Aerosol
- Filters
- High volume fluid samplers
- Electrostatic precipitators
- Settling plates

# Aerosol Sampling

- Efficacy affected by
  - Collection media
  - Relative humidity
  - Desiccation

# BioWatch Program





# Background

## ➤ Rising concerns

- Terrorists could aerosolize biological agents in the air
- Potentially causing thousands of casualties
- Without early detection, the magnitude of the problem might only be revealed as people arrived at hospitals with symptoms

## ➤ Early-warning system

- Department of Homeland Security
- Tested at the winter Olympic, Salt Lake City in 2002
- BioWatcat Program in 2003

# Agencies & Responsibilities

- Department of Homeland Security (DHS)
- Environmental Protection Agency (EPA)
- Center for Disease Control and Prevention (CDC)
- Department of Energy (DOE)
- Department of Defense (DOD)

## DHS

Funding,  
Oversight,  
Coordination

## DOE and DOD

Site selection,  
Sampler unit,  
Tech support

## EPA

Field coordination,  
Site maintenance,  
Sampling operation

## CDC

Sample analysis &  
data management,  
Consequence management  
strategy

# Sampling


- Outdoor
  - Pre-existing EPA air quality station
  - Manually collect samples 24-hr intervals, 7 days/week, 365 days/year
    - Potable sampling unit (PSU)
    - Dry filter unit (DFU)
    - Distributed sampling unit
  - Deliver samples daily to CDC-operated facilities
- Indoor (enhancing & developing the protocols)
  - Suggested methods
    - HEPA socks?
    - Swabs for areas where water pools?
    - Swipes for surfaces, e.g. computer screens, air intake vents?

# Analysis & Responses

- CDC-operated BioWatch facility
  - Federal, State, or local public health laboratories
  - Primary analysis
    - Polymerase chain reaction (PCR)
    - If a biological agent is detected, perform 2<sup>nd</sup>ary test
  - Secondary analysis
    - PCR-verified positive
      - existence of a biological material in the air due to either an intentional or a natural
- Notification by local level with health authorities
- Low enforcement by Federal Bureau Investigation (FBI)

# BioWatch Program Model

As shown below, the BioWatch program has three components: sampling (A), analysis (B), and response (C).



# Evaluations of the Program



**The Mercury News**  
MercuryNews.com

Current: 66°  
67° / 59°

[Complete Forecast](#)

News ☐ Archives ☐ Web ☐ for    
[Back to Home >](#)

Welcome Guest  
[Sign Up](#) | [Sign In](#) | [Member Benefits](#)  
Wednesday, Nov 09, 2005

## News

Posted on Wed, Nov. 02, 2005

### Terror germ detected in Santa Clara County, but threat downplayed

SUBSEQUENT AIR TESTS NEGATIVE  
By Chuck Carroll  
Mercury News

Public health officials are downplaying the threat of a bioterror attack in Santa Clara County, saying it probably came from a natural source.

On Sunday morning, an air monitor detected the presence of the tularemia bacteria, but a test of Public Health. No other air monitors detected the bacteria.

Dr. Marty Fenstersheib, San Jose's health officer, said that the source was likely a natural one, but he did not say which cities were most at risk.

"All of our subsequent tests were negative," he said.

The 31 cities from which the bacteria were detected are listed below.


### Test Results Cited in Delay of Mall Alert CDC Explains Why Local Officials Weren't Told for Days About Bacterium Detection

By Susan Levine and Sari Horwitz  
Washington Post Staff Writers  
Wednesday, October 5, 2005; Page B01

Area health officials were not notified for five days that a potentially dangerous bacterium had been detected because subsequent tests were negative, a CDC official said yesterday.

- EPA concerns
  - Quality assurance and quality control (QAQC)
  - Protocols and response systems

# Detection Goals



# Three Dimensions of Detection

## Viruses

Rotavirus  
Hepatitis A  
Norwalk virus

## Protozoa

Cryptosporidium  
Giardia

## EPA: Indicator organisms

*E. coli*  
Fecal enterococci  
Bacteroides spp.  
Bifidobacteirum spp.  
Clostridium spp.  
Eubacterium spp.  
Background Organism (List of phylogentic groups)

## CDC Category B- Water

*Yersinia enterocolitica*  
*Campylobacter*  
*Listeria monocytogenes*  
Salmonella (Newport and DT104)  
Pathogenic vibrios

Shigella  
Diarrheagenic *E. coli*

**Breadth: Species**

**Function:** antibiotic resistance, virulence

*E. coli* K12

*E. coli* O157:H7

*E. coli* O26:H11

*E. coli* O111:H8

*E. coli* O121:H19

**Depth: Strain**



# Polymerase Chain Reaction (PCR)


- Developed in 1985 by Kary Mullis
- Dr. Mullis received the Nobel Prize in Chemistry in 1993
- PCR is considered as one of the most important discoveries in molecular biology

# What is PCR?

- Enzymatic reaction that makes many copies of DNA from single molecule
- $2^n$  copies of DNA from single molecule where  $n$  = No. of cycles
- So, 35 cycles of PCR would yield  $2^{35}$  copies of DNA


# PCR : Polymerase Chain Reaction

30 - 40 cycles of 3 steps :



**Step 1 : denaturation**


1 minut 94 °C



**Step 2 : annealing**

45 seconds 54 °C

forward and reverse  
primers !!!



**Step 3 : extension**


2 minutes 72 °C

only dNTP's

(Andy Vierstraete 1999)

# Applications

- Allows for amplification and identification of specific gene sequences from specific organisms (i.e. bacteria, parasites, viruses, etc.) that would otherwise be undetectable
- Provides information for HAZ ID
- Provides information for sources (HAZ)
- Used for Forensic studies (Bioterrorism)
- Used for molecular epidemiology (outbreak investigation)




Source

Timothy D. Read *et al.* Comparative genome sequencing for discovery of novel polymorphisms in *Bacillus anthracis*. *Science* 296:2028-2033 (2002).



# Esp marker



- Host specific genetic marker  
Differentiate Human/non human pollution
- Present in Indicator bacteria (Enterococci)
- Traditional MF method and PCR used to detect the gene
- A culture dependent library independent method

# Methods

**Membrane filtration**



**Incubation (24hrs)**



**Method 1  
conventional**

**Method 2  
Rapid**

**Method 3  
Rapid**



**Enrichment (2 hrs)**

**Enrichment (2hrs)**

**Cell concentration  
(30 minutes)**



**DNA extraction  
(3.5hrs)**

**Cell concentration  
(30 minutes)**




**PCR (3hrs)**

# Membrane Filtration (MF)





# Enterococci (After Incubation)



**Colonies that  
have a blue  
halo are  
counted as  
enterococci.**

# Enrichment





# Polymerase Chain Reaction (PCR)

- Enzymatic reaction that makes many copies of DNA from single molecule
- Forward primer, which is specific for the *E. faecium esp* gene is: (5'-TAT GAA AGC AAC AGC ACA AGT T-3')
- Conserved reverse primer (5'-ACG TCG AAA GTT CGA TTT CC-3')

# Gel loading



# Esp detection after PCR



qReal Time PCR: What can it do for my project?

## Real time PCR : in situ detection

Detect PCR product fluorescently in each well of PCR plate.


Fast PCR screening without gels.

Quantify amount of PCR product at each cycle.

Detect presence or quantify fraction of sample made up by particular species using species specific primers.

Use threshold detection for relative abundance.


## REAL TIME PCR using Taqman probe method






## Polymerization


(R) = Reporter  
(Q) = Quencher



## Strand displacement



## Cleavage



## Polymerization completed




**Well-designed TaqMan probes require very little optimization.**

**In addition, they can be used for multiplex assays by designing each probe with a spectrally unique fluor/quench pair.**

**However, TaqMan probes can be expensive to synthesize, with a separate probe needed for each DNA target being analyzed.**

**Each probe is ~\$ 400.00**

# Real-time QPCR using molecular beacon



# Factors Important in Developing Detection Tools

## Factors to Consider in Developing Parallel Detection Tools for Microorganisms

### TARGET SELECTION

#### Bacteria

*Bacillus anthracis*  
*E. coli* O157:H7  
Etc.

#### Protozoa

*Cryptosporidium*  
*Giardia*  
Etc.

#### Viruses

Norwalk virus  
SARS  
Etc.

### PROBE DESIGN

#### Gene Target

Small subunit rRNA  
Functional genes  
Intragenic spacer region  
Virulence factors

#### Specificity

Probe size  
Perfect match/mismatch strategy

#### Resolution

Genus & species level  
Strain fingerprinting

### MICROARRAY PLATFORM

#### Glass slide arrays

In-house printed arrays  
Commercial e.g., Agilent

#### In situ synthesized arrays

e.g., Affymetrix  
NimbleGen  
Xeotron

#### Gel-pad arrays

#### Diagnostic arrays

e.g., Motorola's eSensor  
Lab-on-a-chip systems

### SAMPLE PROCESSING

#### Matrices

Air, water, food  
Body fluid, tissue, soil  
Others

#### Sample Concentration

Micro- & Ultrafiltration  
Centrifugation, others

#### Target Enrichment & Amplification

Bead-based separation  
Polymerase chain reaction  
Biological amplification on media

### VALIDATION & DATA ANALYSIS

#### Detection

Identification  
Reliability  
False positives  
False negatives  
Detection limit  
Signal amplification

#### Quantification

Dissociation curves  
Standard mixtures  
Chip-based PCR

### USABILITY

#### Upgradability

Database growth


#### User-friendliness

Simplicity  
Speed  
Portability  
Robustness

#### Competitiveness








Cost  
Niche  
Competing technologies

# Limit of Detection



Cells or copies per reaction well (volume = a few  $\mu\text{l}$  to 1 mL)

# Comparison of 6 PCR-based and 1 Antibody-based Assay

	PCR Chip Hand-Held	Arrow-Tech Bio-Seq (HANAA)	Idaho RAZOR	PCR-Chip Bench top	Idaho R.A.P.I.D.	BioTrove NTCycler	BioVeris M1M analyzer
Base Unit	\$1,000	\$31,000	\$35,000	\$8,000	\$55,000	\$100,000	\$70,000
Consumable Price (U)	\$10	\$38	\$18	\$100		\$400	\$1,000
Time	30 min	20-30 min	30 min	30 min	30 min	3.5 hr	90 min
# of Samples	14	6	12	56	10	48	96
# of Targets	50	1-4	12	50	10	18-56	1
Sample volume	10-500 nl		5-20 µl	10-500 nl	5-20 µl	33 nl	100 µl
Quantitation	Yes	No	No	No	Yes	Yes	No
Sensitivity	5 copies	5 copies	5 copies	5 copies	5 copies	5 copies	1000 cells
Portability (Weight)	< 1 lb	6.5 lb	9 lb	30 lb	50 lb	100 lb	100 lb
	<b>Goal*</b> 						

**Thank You**